

LifeGel Plates preparation for cell-based assays.

Materials required:

- 48-, 96- or 384-well LifeGel plate [RRLG96; RRLG48; RRLG384]
- Pipette/multichannel pipettor
- Complete culture medium



1. Equilibration

Remove LifeGel plates from the refrigerator 24 or 48 hours before seeding. Allow them to equilibrate at room temperature for approximately 1 hour.

2. Gel conditioning


Aspirate the PBS from each well*. Add 400 μ l of culture medium per well for 48-well LifeGel plates, 150 μ l for 96-well LifeGel plates or 60 μ l for 384-well plates to condition the gel. Incubate the plates at 37°C in a humidified incubator.

3. Cell seeding

After 24 or 48 hours of conditioning, remove the plates from the incubator. Aspirate the conditioning medium from the wells*.

Add 400 μ l (48-well), 150 μ l (96-well) or 60 μ l (384-well) of fresh culture medium containing cell suspension to each well.

Return the seeded LifeGel plates to the 37°C incubator for cell culture.



* Be careful when aspirating PBS or conditioning medium to avoid disturbing the gel.