

IMMUNOFLUORESCENCE STAINING OF ACTIN, TUBULIN AND NUCLEI USING PANC-1 AND HEK293 CELLS GROWN ON LIFE GEL



IMMUNOFLUORESCENCE STAINING OF ACTIN, TUBULIN AND NUCLEI USING PANC-1 AND HEK293 CELLS GROWN ON LIFEHEL

PROCEDURE

LifeGel on 15 well μ -slide from IBIDI

This procedure was optimized for cells growing on LifeGel on 15 well μ -slide from IBIDI.

1. Fix cells using 4% formaldehyde (10% formalin solution) – 40 μ L per well and incubate for 15 minutes at room temperature.
2. Wash cells 3x with 40 μ L of PBS.
3. Perforate cells using 0,1% Triton X-100 in PBS with 10% FBS. Use 40 μ L of solution per well and incubate for 15 minutes at room temperature.
4. Wash cells twice with PBS.
5. Change PBS for 40 μ L of working solution of primary antibody (1:1000 dilution) and incubate plate overnight at 4°C.
6. Next day wash cells 3 times with PBS (wait 5 minutes between each washing).
7. Prepare mix of working solution containing secondary antibody (1:500), phalloidin (1:500) and Hoechst (10 μ g/mL) in PBS. Add 30 μ L per well and incubate plate for 3 hours at room temperature in darkness.
8. Wash cells with PBS, 3 times (wait 5 minutes between each washing).
9. Then take images of brightfield and fluorescence using fluorescence microscope and objectives characterized by working distance not less than 1,5 mm – objectives with higher magnifications like 10x and 20x are preferable for cytoskeleton training.

LifeGel on 96 well plate

1. Fix cells using 4% formaldehyde (10% formalin solution) – 100 μ L per well and incubate for 15-30 minutes at room temperature (the bigger 3D structures, the longer time of incubation is recommended).
2. Wash cells 3x with 150 μ L of PBS.
3. Perforate cells using 0,1% Triton X-100 in PBS with 10% FBS. Use 100 μ L of solution per well and incubate for 15-30 minutes at room temperature (time similar as for cell fixation).
4. Wash cells twice with PBS.
5. Change PBS for 100 μ L of working solution of primary antibody (1:500 dilution) and incubate plate overnight at 4°C.
6. Next day wash cells 3 times with PBS (wait 5 minutes between each washing).
7. Prepare mix of working solution containing secondary antibody (1:250), phalloidin (1:500) and Hoechst (10 μ g/mL) in PBS. Add 100 μ L per well and incubate plate for 3 hours at room temperature in darkness.
8. Wash cells with PBS, 3 times (wait 5 minutes between each washing).
9. Then take images of brightfield and fluorescence using fluorescence microscope and objectives characterized by working distance not less than 2 mm – objectives with higher magnifications like 10x and 20x are preferable for cytoskeleton imaging.




Fig. 1. PANC-1 spheroids immunostaining on LifeGel. Objective 10x, optovar lenx 2,5x.




Fig. 2. HEK293 cells immunostaining on LifeGel. Objective 10x, optovar lens 2,5x.

Find out more on LifeGel and another application notes at

For research use only.

Real-Research products are intended for research use only, not for diagnostic or therapeutic purpose.

Real-Research sp. z o.o. ul. prof. Michała Bobrzyńskiego 14, 30-348 Kraków
+48 603-984-577 | sale@real-research.com | www.real-research.com

